

Ez SHOT ID

HIL

Easy to Use Auto Focus HD intraoral camera

10

Ju Sto Ex Stron Inte

CE

Intraoral Camera

-IC-WHT60 -IC-WHW61 -IC-EZAU81 Curing Light

Apex Locator

Easy to Use Auto Focus HD intraoral camera

1 Features of Ez SHOT

Auto focus

HD	

HD image Can get high definition image without distortion.



Can get clear image with auto focus function in oral cavity.

Adjustable LED brightness To avoid image quality loss by light reflection, can adjust LED brightness.



Black & White function Can take and send analog X-ray film image to PC.



No need additional driver to install the device.

Easy to capture Easy to image capture by camera only without extra device like foot pedal.


Software compatibility Can be use camera with most of dental imaging software. ^t It would need twain driver for some of dental imaging software.

2 Specification

Model name		IC-EZAU81
Camera	Operation voltage	5V
	Size	Approx. 202 x 27 x 29 mm
	Weight	Approx. 55g
	Resolution	1.3 Mega pixel
	Sensor/Lens type	CMOS / Liquid lens
USB cable	Length	Approx. 3M
Disposable cover	Size	Approx. 32 x 100 mm
	Quantity	Approx. 100 pcs by weight x 2ea
	Raw material	Polyethylene

3 Minimum System Requirement

PC system:	2.0 Ghz CPU, 1GB RAM	
OS:	Windows XP service pack 3	
USB port:	USB 2.0	



Upper images are take by real product

4 DrsViewer2

Image Viewer included.





Good Doctors Co., Ltd.

#208, B-dong, 283 Bupyeong-daero, Bupyeong-gu (Woolim Lions Valley, Cheongchoen-dong), Incheon, 21315, Republic of Korea
Tel: +82 32 424 6325 Fax: +82 32 424 6326
E-mail: info@gooddrs.com Website: www.gooddrs.com

Gwangju Branch

#203, Sauphwa 2–dong, 333 Cheomdangwagi–ro, Buk–gu (Gwangju Techno Park, Daechon–dong), Gwangju, 61008, Republic of Korea Tel: +82 62 602 7928 Fax: +82 62 602 7426

